

# Gracey Curettes

Nordent's Gracey curettes are designed just like Dr. Gracey intended them to be. Using Nordent Gracey curettes will reduce tissue trauma, minimize lateral pressure and reduce the number of strokes required during a procedure.



**All Nordent Gracey patterns are available in:**

- Standard
- Rigid ("R" suffix)
- Long ("L" suffix)
- Mini-Long ("ML" suffix)

## Gracey Designs

**Standard**

Standard terminal shank and blade.

**Rigid**



"Rigid" terminal shank with a slightly longer and wider blade than a Standard Gracey. A good choice for periodontal surgery.

**Long**

The terminal shank is 3mm longer than the Standard Gracey with a Standard blade. For improved access in deep periodontal pockets.

**Mini-Long**

The terminal shank is 3mm longer than the Standard Gracey with a shorter and narrower "Mini" blade. Excellent for periodontal maintenance in deep pockets



## Gracey Application Chart

Use this color-coded chart to select the correct Gracey patterns for your application.

**Anteriors**

-  All Surfaces  
Gracey Patterns 1-2, 3-4 or 5-6

**Premolars and Molars**

-  Buccal & Lingual Surfaces  
Gracey Patterns 7-8 or 9-10
-  Mesial Surfaces  
Gracey Patterns 11-12 or 15-16
-  Distal Surfaces  
Gracey Patterns 13-14 or 17-18

# Gracey Curettes

Nordent Standard Gracey curettes feature a blade face curve with a true radius that conforms to tooth anatomy better than any other brand. This means that Nordent Standard Gracey curettes are easier to adapt, require fewer strokes, and reduce tissue trauma better than any other brand. Nordent Standard Gracey curettes have a 1.1mm terminal shank diameter and blades that are 0.9mm wide and 4mm long.

## Standard



**Gracey #1-2** Short contra-angle for incisors and cuspids. Terminal shank reach is 11mm.

Handle Selection

CESCGR1-2 shown  RESCGR1-2  ESCGR1-2  RSCGR1-2  SCGR1-2



**Gracey #3-4** Short contra-angle for incisors and cuspids. Terminal shank reach is 11mm.

Handle Selection

CESCGR3-4 shown  RESCGR3-4  ESCGR3-4  RSCGR3-4  SCGR3-4



**Gracey #5-6** Medium shank and slight contra-angle for bicuspid and molars. Terminal shank reach is 14mm.

Handle Selection

CESCGR5-6 shown  RESCGR5-6  ESCGR5-6  RSCGR5-6  SCGR5-6



**Gracey #7-8** Medium contra-angle for facial and lingual surfaces of bicuspid and molars. Terminal shank reach is 12mm.

Handle Selection

CESCGR7-8 shown  RESCGR7-8  ESCGR7-8  RSCGR7-8  SCGR7-8



**Gracey #9-10** Long contra-angle for molars and less accessible root surfaces. Terminal shank reach is 12mm.

Handle Selection

CESCGR9-10 shown  RESCGR9-10  ESCGR9-10  RSCGR9-10  SCGR9-10

## Nordent Tip

Gracey curettes are also known as "Area Specific Curettes". All Gracey's have one cutting edge on the lower side of the blade and a rounded toe. When re-sharpening a Gracey blade sharpen the lower edge and the toe only. Visit [www.nordent.com/instrenew](http://www.nordent.com/instrenew) to see a video demonstration of a simple and effective way sharpen any brand and any pattern of scalers and curettes.

A complete set of Gracey curettes would include one for Anteriors, one for the buccal and lingual surfaces of pre-molars and molars, one for the mesial surfaces of pre-molars and molars and one for the distal surfaces of premolars and molars.

# Gracey Curettes

Standard

HYGIENE



**Gracey #11-12** All mesial surfaces of bicuspid and molars. The terminal shank reach is 8mm.

Handle Selection  CESCGR11-12 shown  RESCGR11-12  ESCGR11-12  RSCGR11-12  SCGR11-12



**Gracey #13-14** All distal surfaces of bicuspid and molars. The terminal shank reach is 9mm.

Handle Selection  CESCGR13-14 shown  RESCGR13-14  ESCGR13-14  RSCGR13-14  SCGR13-14



**Gracey #11-14** Mesial Distal Combination combines the mesial tip of the Gracey 11 and the distal Gracey 14 tip to allow completion of the facial or lingual sextant without changing instruments.

Handle Selection  CESCGR11-14 shown  RESCGR11-14  ESCGR11-14  RSCGR11-14  SCGR11-14



**Gracey #12-13** Mesial Distal Combination combines the mesial tip of the Gracey 12 and the distal Gracey 13 tip to allow completion of the facial or lingual sextant without changing instruments.

Handle Selection  CESCGR12-13 shown  RESCGR12-13  ESCGR12-13  RSCGR12-13  SCGR12-13



**Gracey #15-16** Modified version of the Gracey 11-12 allows better access and a more comfortable hand and finger position for all mesial surfaces of bicuspid and molars. The terminal shank reach is 9mm.

Handle Selection  CESCGR15-16 shown  RESCGR15-16  ESCGR15-16  RSCGR15-16  SCGR15-16



**Gracey #17-18** Modified version of the Gracey 13-14 allows better access and a more comfortable hand and finger position for all distal surfaces of bicuspid and molars. The terminal shank reach is 7mm.

Handle Selection  CESCGR17-18 shown  RESCGR17-18  ESCGR17-18  RSCGR17-18  SCGR17-18

# Gracey Curettes

## Long Reach

Nordent Long Reach Gracey curettes have an extended terminal shank that is 3mm longer than standard. They feature our standard blade face curve with a true radius that conforms to tooth anatomy better than any other brand. This means that Nordent Standard Gracey curettes are easier to adapt, require fewer strokes, and reduce tissue trauma better than any other brand. Nordent Standard Gracey curettes have a 1.1mm terminal shank diameter and blades that are 0.9mm wide and 4mm long.



### Long Reach Gracey #1-2 Short contra-angle for incisors and cuspids.

Handle Selection  CESCGR1-2L shown  RESCGR1-2L  ESCGR1-2L  RSCGR1-2L  SCGR1-2L



### Long Reach Gracey #3-4 Short contra-angle for incisors and cuspids.

Handle Selection  CESCGR3-4L shown  RESCGR3-4L  ESCGR3-4L  RSCGR3-4L  CESCGR3-4L



### Long Reach Gracey #5-6 Medium shank and slight contra-angle for bicuspid and molars.

Handle Selection  CESCGR5-6L shown  RESCGR5-6L  ESCGR5-6L  RSCGR5-6L  SCGR5-6L



### Long Reach Gracey #7-8 Medium contra-angle for facial and lingual surfaces of bicuspid and molars.

Handle Selection  CESCGR7-8L shown  RESCGR7-8L  ESCGR7-8L  RSCGR7-8L  SCGR7-8L



### Long Reach Gracey #11-12 All mesial surfaces of bicuspid and molars.

Handle Selection  CESCGR11-12L shown  RESCGR11-12L  ESCGR11-12L  RSCGR11-12L  SCGR11-12L



### Long Reach Gracey #13-14 All distal surfaces of bicuspid and molars.

Handle Selection  CESCGR13-14L shown  RESCGR13-14L  ESCGR13-14L  RSCGR13-14L  SCGR13-14L



### Long Reach Gracey #15-16 Modified version of the Long Reach Gracey 11-12 allows better access and a more comfortable hand and finger position for all mesial surfaces of bicuspid and molars.

Handle Selection  CESCGR15-16L shown  RESCGR15-16L  ESCGR15-16L  RSCGR15-16L  SCGR15-16L



### Long Reach Gracey #17-18 Modified version of the Long Reach Gracey 13-14 allows better access and a more comfortable hand and finger position for all distal surfaces of bicuspid and molars.

Handle Selection  CESCGR17-18L shown  RESCGR17-18L  ESCGR17-18L  RSCGR17-18L  SCGR17-18L

# Gracey Curettes

## Mini Blade - Long Reach

Nordent Mini Blade - Long Reach Gracey curettes have a shorter, narrower blade and an extended terminal shank that is 3mm longer than standard. They feature our blade face curve with a true radius that conforms to tooth anatomy better than any other brand. This means that Nordent Gracey curettes are easier to adapt, require fewer strokes, and reduce tissue trauma better than any other brand. Nordent Mini Blade - Long Reach Gracey curettes have a 1mm terminal shank diameter and blades that are 0.7mm wide and 3mm long.



**Mini Blade - Long Reach Gracey #1-2** Short contra-angle for incisors and cuspids.

Handle Selection  CESCGR1-2ML shown  RESCGR1-2ML  ESCGR1-2ML  RSCGR1-2ML  SCGR1-2ML



**Mini Blade - Long Reach Gracey #3-4** Short contra-angle for incisors and cuspids.

Handle Selection  CESCGR3-4ML shown  RESCGR3-4ML  ESCGR3-4ML  RSCGR3-4ML  SCGR3-4ML



**Mini Blade - Long Reach Gracey #5-6** Medium shank and slight contra-angle for bicuspid and molars.

Handle Selection  CESCGR5-6ML shown  RESCGR5-6ML  ESCGR5-6ML  RSCGR5-6ML  SCGR5-6ML



**Mini Blade - Long Reach Gracey #7-8** Medium contra-angle for facial and lingual surfaces of bicuspid and molars.

Handle Selection  CESCGR7-8ML shown  RESCGR7-8ML  ESCGR7-8ML  RSCGR7-8ML  SCGR7-8ML



**Mini Blade - Long Reach Gracey #11-12** All mesial surfaces of bicuspid and molars.

Handle Selection  CESCGR11-12ML shown  RESCGR11-12ML  ESCGR11-12ML  RSCGR11-12ML  SCGR11-12ML



**Mini Blade - Long Reach Gracey #13-14** All distal surfaces of bicuspid and molars.

Handle Selection  CESCGR13-14ML shown  RESCGR13-14ML  ESCGR13-14ML  RSCGR13-14ML  SCGR13-14ML



**Mini Blade - Long Reach Gracey #15-16** Modified version of the Mini Blade - Long Reach Gracey 11-12 allows better access and a more comfortable hand and finger position for all mesial surfaces of bicuspid and molars.

Handle Selection  CESCGR15-16ML shown  RESCGR15-16ML  ESCGR15-16ML  RSCGR15-16ML  SCGR15-16ML



**Mini Blade - Long Reach Gracey #17-18** Modified version of the Mini Blade - Long Reach Gracey 13-14 allows better access and a more comfortable hand and finger position for all distal surfaces of bicuspid and molars.

Handle Selection  CESCGR17-18ML shown  RESCGR17-18ML  ESCGR17-18ML  RSCGR17-18ML  SCGR17-18ML

# Gracey Curettes

## Rigid

Nordent Rigid Gracey curettes feature a Rigid shank and a longer reach blade and are well suited to periodontal and surgical application. Nordent Rigid Gracey curettes have a 1.3mm terminal shank diameter and blades that are 0.9mm wide and 6mm long.



**Rigid Gracey #1-2** Short contra-angle for incisors and cuspids. Terminal shank reach is 11mm.

Handle Selection

CESC1R-2R shown  RESC1R-2R  ESC1R-2R  RSC1R-2R  SC1R-2R



**Rigid Gracey #3-4** Short contra-angle for incisors and cuspids. Terminal shank reach is 11mm.

Handle Selection

CESC3R-4R shown  RESC3R-4R  ESC3R-4R  RSC3R-4R  SC3R-4R



**Rigid Gracey #5-6** Medium shank and slight contra-angle for bicusps and molars. Terminal shank reach is 14mm.

Handle Selection

CESC5R-6R shown  RESC5R-6R  ESC5R-6R  RSC5R-6R  SC5R-6R



**Rigid Gracey #7-8** Medium contra-angle for facial and lingual surfaces of bicusps and molars. Terminal shank reach is 12mm.

Handle Selection

CESC7R-8R shown  RESC7R-8R  ESC7R-8R  RSC7R-8R  SC7R-8R



**Rigid Gracey #9-10** Long contra-angle for molars and less accessible root surfaces. Terminal shank reach is 12mm.

Handle Selection

CESC9R-10R shown  RESC9R-10R  ESC9R-10R  RSC9R-10R  SC9R-10R



**Rigid Gracey #11-12** All mesial surfaces of bicusps and molars. The terminal shank reach is 8mm.

Handle Selection

CESC11R-12R shown  RESC11R-12R  ESC11R-12R  RSC11R-12R  SC11R-12R



**Rigid Gracey #13-14** All distal surfaces of bicusps and molars. The terminal shank reach is 7mm.

Handle Selection

CESC13R-14R shown  RESC13R-14R  ESC13R-14R  RSC13R-14R  SC13R-14R