

Explorers

Explorers are used to examine teeth for decay (caries), calculus and other abnormalities. They have thin, flexible tips that taper to a sharp point and come in various single end and double end patterns. All Nordent Explorer tips are made from spring tempered stainless steel and are machine ground to exacting tolerances. Each tip is then hand sharpened and formed into its final shape by expert craftsmen. With Nordent, you can be assured of the highest level of quality, consistency, and the sharpest tips in the industry.

Single End



#23

Handle Selection CEEX23 *shown* REEX23 EX23



#408

Handle Selection CEEXS408 *shown* REEXS408 EXS408



#3A

Handle Selection CEEX3A *shown* REEX3A EX3A



#6

Handle Selection CEEX6SE *shown* REEX6SE EX6SE



#6A

Handle Selection CEEX6A *shown* REEX6A EX6A



#6XL

Handle Selection CEEX6XL *shown* REEX6XL EX6XL



#9

Handle Selection CEEXS9 *shown* REEXS9 EXS9



#17

Handle Selection CEEX17 *shown* REEX17 EX17



#23-6

Handle Selection CEEX23-6 shown REEX23-6 EX23-6



#5

Handle Selection CEEX5 shown REEX5 EX5



Tufts #17-23

Handle Selection CEEXU17-23 shown REEXU17-23 EXU17-23



#6DE

Handle Selection CEEX6DE shown REEX6DE EX6DE

Double End - Endodontic



#DG16

Handle Selection CEEX16 shown REEX16 EX16



#DG16-23

Handle Selection CEEX16-23 shown REEX16-23 EX16-23

Explorers

Double End



#2 Pig Tail

Handle Selection CEEX2 shown REEX2 EX2



#2A

Handle Selection CEEX2A shown REEX2A EX2A



#2XL

Handle Selection CEEX2XL shown REEX2XL EX2XL



#3DE

Handle Selection CEEX3DE shown REEX3DE EX3DE



#3 Cow Horn

Handle Selection CEEX3CH shown REEX3CH EX3CH



ODU - Old Dominion University #11-12

Handle Selection CEEX11-12 shown REEX11-12 EX11-12



ODU - Old Dominion University #11-12L With extended terminal shank.

Handle Selection CEEX11-12L shown REEX11-12L EX11-12L