

Periodontal Probes

Probes are used to measure periodontal pocket depth in millimeter increments. All Nordent Probe tips are made from spring tempered stainless steel and machine ground to exacting tolerances. The millimeter measurements are precisely cut into the Probe tips. Color coding is applied to Probes under extreme temperatures so that the color penetrates the surface of the tip. Probe tips undergo additional tight tolerance grinding to achieve the desired readability for each marking. All tips are then hand formed into their final shape by expert craftsmen.

Single End



N12 Color Coded University of North Carolina Screening Probe 1-2-3-4-5-6-7-8-9-10-11-12 mm markings.

Handle Selection CEPCN12 *shown* REPCN12 PCN12



N15 Color Coded University of North Carolina Screening Probe 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15 mm markings.

Handle Selection CEPCN15 *shown* REPCN15 PCN15



WHO Color Coded "PSR" Ball Tip Screening Probe 3.5-5.5-8.5-11.5 mm markings.

Handle Selection CEPCWHO *shown* REPCWHO PCWHO



N22 Color Coded Probe 2-4-6-8-10-12 mm markings.

Handle Selection CEPCN22 *shown* REPCN22 PCN22



N33 Color Coded Probe 3-6-9-12 mm markings.

Handle Selection CEPCN33 *shown* REPCN33 PCN33



N8-11 Color Coded Probe 3-6-8-11 mm markings.

Handle Selection CEPCN8-11 *shown* REPCN8-11 PCN8-11

Periodontal Probes

Single End



W Non-Color Coded Williams Probe 1-2-3-5-7-8-9-10 mm markings.

Handle Selection CEPW shown REPW PW



OW Color Coded Williams Probe 1-2-3-5-7-8-9-10 mm markings.

Handle Selection CEPCOW shown REPCOW PCOW

Double End



GF-W Non-Color Coded Goldman Fox (flat) - Williams (round) Probe 1-2-3-5-7-8-9-10 mm markings both tips.

Handle Selection CEPGF-W shown REPGF-W PGF-W



N116 Non-Color Coded Offset Probe 1-2-3-5-7-8-9-10 mm markings. "Mirror Image" shanks are curved 25° to allow access to posterior areas.

Handle Selection CEPN116 shown REPN116 PN116



N116 Color Coded Offset Probe 1-2-3-5-7-8-9-10 mm markings. "Mirror Image" shanks are curved 25° to allow access to posterior areas.

Handle Selection CEPCN116 shown REPCN116 PCN116



N136 Color Coded Offset Probe 3-6-9-12 mm markings. "Mirror Image" shanks are curved 25° to allow access to posterior areas.

Handle Selection CEPCN136 shown REPCN136 PCN136

Periodontal Probes

Double End -Bifurcation



NAB2 Color Coded Nabors Furcation Probe 3-6-9-12 mm markings .

Handle Selection



CEPCNAB2 shown



REPCNAB2



PCNAB2



NAB2 Non-Color Coded Nabors Furcation Probe no markings.

Handle Selection



CEPNAB2 shown



REPAB2



PNAB2