

# Handles















## DURALite® ColorRings™ Handle



Our newest and most advanced ergonomic handle design! DuraLite ColorRings has the same large 3/8" (9.5 mm) diameter as our popular DuraLite Round design. The handle design has a more gradual taper from the tip to the grip so that your fingers are even more comfortable closer to the tip. In addition, we've extended the ControlRing grip by 50% to give you maximum slip-free gripping power over a larger portion of the handle. The handle has a gentle texture to reduce glare and is 100% stainless steel to make it easy to clean and sterilize by any method. The DuraLite ColorRing handle incorporates four color ring positions that are below the finger grip surface allowing you to better identify, organize and customize your instrumentation. Our innovative design makes customizing the DuraLite ColorRing handle simple and easy. With 14 Colors available...How will you color your world?

## ColorRings™

ColorRings come in 14 different colors and are available in single color packages of 48-rings. Use the item codes below to select your colors.

													
RING-BL Blue	RING-TE Teal	RING-FS Fuchsia	RING-YW Yellow	RING-LG Lt Green	RING-WI White	RING-RD Red	RING-GN Green	RING-OR Orange	RING-PR Purple	RING-BR Brown	RING-BK Black	RING-GR Gray	RING-PK Pink

## DURALite® Round Handle



Our DuraLite Round has a large 3/8" (9.5 mm) diameter and is extremely lightweight for maximum tactile sensitivity. Our exclusive ControlRing design gives you the ultimate in grip and control. It requires less finger pressure resulting in less hand and finger fatigue. The handle has a gentle texture to reduce glare and is 100% stainless steel to make it easy to clean and sterilize by any method.

## DURALite® HEXagonal



Our original ERGONOMIC handle that changed the industry! Duralite HEX has a large 7/16" (11 mm) diameter that comfortably positions your hand and fingers in the "Modified Pen Grasp" position. The handle has a gentle texture to reduce glare and is 100% stainless steel to make it easy to clean and sterilize by any method.

## Medium Round



The Medium Round #4 handle has a moderate 5/16" (8mm) diameter with a segmented knurled grip to provide excellent control. The handle is 100% stainless steel and the finish is bright electro-polished to make it easy to clean and sterilize by any method.

## Standard



This is Nordent's original solid handle design. It has a 1/4" (6.4 mm) diameter with a segmented knurled grip to provide excellent control with this slimmer handle design. The handle is 100% stainless steel and the finish is bright electro-polished to make it easy to clean and sterilize by any method.