

# Single End Scalers

Single end scalers are an excellent choice for pedodontic practices to remove deposits from supragingival surfaces and in orthodontic practices to remove cement. Single end scalers should be used when a double end instrument is not suitable or required. Nordent will produce any of our hygiene tip patterns in a single end configuration upon request.



**Sickle #5** A moderate size, medium-reach sickle with a blade width of 0.9mm that tapers to a point.

Handle Selection  CESC#5 *shown*  RESC#5  SC#5



**Sickle #109** A moderate size, long-reach sickle with a blade width of 0.9mm that tapers to a point.

Handle Selection

CESC109  
 RESC109  
 SC109



**Towner U15** A large size, long-reach sickle with a blade width of 1mm that tapers to a point.

Handle Selection

CESC#15  
 RESC#15  
 SC#15



**Ivory CI-1** A moderate size, medium-reach sickle with a blade width of 2mm that tapers to a point.

Handle Selection

CESC#1  
 RESC#1  
 SC#1



**Younger-Good #15** A small size, short-reach sickle on a straight shank with a blade width of 1mm that tapers to a point.

Handle Selection

CESC#15  
 RESC#15  
 SC#15



**Univ. of So. Cal. #128** A moderate size, medium-reach sickle on a straight shank with a blade width of 1.2mm that tapers to a point.

Handle Selection

CESC#128  
 RESC#128  
 SC#128



**Whiteside #2** A small size, short-reach sickle on a straight shank with a 1.5mm blade that tapers to a point.

Handle Selection

CESC#2  
 RESC#2  
 SC#2



**Jacquette J1S** A small size, straight blade Jacquette with a 0.9mm blade width that tapers to a point. This tip is also known as the Jacquette #33.

Handle Selection

CESC#1S  
 RESC#1S  
 SC#1S



**"B" Scaler** A moderate size, straight blade with a 0.5mm width that has two sharp edges on the face and the back of the blade. This enables the blade to be used in a push or pull stroke.

Handle Selection

CESC#B  
 RESC#B  
 SC#B