

# Goldman-Fox

The Goldman-Fox series of instruments includes universal curettes, Jacquettes, sickles and hoes.

## Scalers & Curettes



**GF #1 Anterior Jacquette** The blades are 5mm long and taper to a point from a 1mm width.

Handle Selection  CESC GF1 shown  RESCGF1  ESCGF1  RSCGF1  SCGF1



**GF #2 Anterior Curette** Designed like a Gracey medium contra-angle with an open blade curve and a rigid shank. The blade is 1.1mm wide with a round toe.

Handle Selection  CESC GF2 shown  RESCGF2  ESCGF2  RSCGF2  SCGF2



**GF #3 Universal Curette** The terminal shank is 1.1mm in diameter and 10mm long. The blade width is 1mm.

Handle Selection  CESC GF3 shown  RESCGF3  ESCGF3  RSCGF3  SCGF3



**GF #4 Universal Curette** The terminal shank is 1.1mm in diameter and 12mm long. The blade width is 1mm.

Handle Selection  CESC GF4 shown  RESCGF4  ESCGF4  RSCGF4  SCGF4



**GF #5 Hoe** Mirror image hoes that are 1.5mm wide. For facial and lingual surfaces of posterior teeth.

Handle Selection  CESC GF5 shown  RESCGF5  ESCGF5  RSCGF5  SCGF5



**GF #6 Hoe** The blades are 1.5mm wide and are well suited for anterior areas.

Handle Selection  CESC GF6 shown  RESCGF6  ESCGF6  RSCGF6  SCGF6



**GF #21 Anterior Scaler** Combines a small Jacquette that has a 3.5mm length and a small sickle that has a 1.2mm blade that tapers to a point.

Handle Selection  CESC GF21 shown  RESCGF21  ESCGF21  RSCGF21  SCGF21