


Implant Maintenance

Introducing ImplaMate™ “The Ultimate in Implant Scaler Technology!” Until now clinicians have relied on a limited selection of instruments with thick plastic tips to clean implants. This is because regular scaler tips made of metals that are dissimilar to titanium (like stainless steel) could scratch the implant and possibly cause rejection. Our new series of ImplaMate scalers work safely with implants because they are made with the same material... solid titanium. Even better, ImplaMate scalers are available in a wide range of familiar patterns that are easy to adapt with no “learning curve”. Produced in our DuraLite® ColorRings™ handle, ImplaMate Implant Scalers will also enable your practice to identify and organize implant hygiene set-ups.



IMPLA-MATE Universal Curette/Langer 1-2

Handle Selection  CEISLN1-2 shown




IMPLA-MATE Universal Curette/Langer 3-4

Handle Selection  CEISLN3-4 shown




IMPLA-MATE Universal Curette/Langer 5-6

Handle Selection  CEISLN5-6 shown




IMPLA-MATE Sickle 6-7

Handle Selection  CEISN67 shown




IMPLA-MATE Posterior Sickle 204S

Handle Selection  CEIS204S shown



IMPLA-MATE Universal Curette/Barnhart 5-6

Handle Selection  CEISBH5-6 shown