

# Orthodontic Instruments



**Band Pusher - Black's #6-7** Two identical oval-shaped, serrated blades set at opposing angles. The blades measure 3.5mm wide and 1.5mm thick.

Handle Selection  CECNB6-7  RECNB6-7  ECNB6-7  CNB6-7



**Band Pusher - Scaler #114**

Handle Selection  EOTN114  OTN114



**Scaler - Ligature Director #115**

Handle Selection  EOTN115  OTN115




**Band Pusher - Scaler #119**

Handle Selection  EOTN119  OTN119



**Ligature Director #120 single end**

Handle Selection  OTN120



**Elastic Placement Hemostat #101X**

The blades have a unique cross cut serration at the tip for a secure grasp. 4 3/4" / 120mm