


Hemostats

Used as a multipurpose instrument in dental practices to grasp and move small objects. Hemostat jaws have interlocking teeth and are designed to clamp off blood vessels during surgical procedures.

Mosquito



Mosquito #101 Straight jaws (4-3/4" / 120mm)

H101



Mosquito #102 Curved jaws (4-3/4" / 120mm)


H102

Kelly



Kelly #103 Straight jaws (5-3/4" / 145mm)

H103



Kelly #104 Curved jaws (5-3/4" / 145mm)

H104