

Luxation Blades & Periostomes

Luxation Blades

Designed to be easily inserted along the root surface and cut the periodontal ligament prior to extraction. Using Nordent Luxation Blades will minimize tissue trauma and promote quicker healing. Features: Long, Slender Profile, Ultra-Sharp blades, Titanium Nitride Coating to retain sharpness longer, Comfortable Form-Fit Stainless Steel Handle.

#3CT
Curved 3mm
blade width.

ELX3CT



#3ST
Straight 3mm blade width.

ELX3ST



#5CT
Curved 5mm blade width.

ELX5CT



#5ST
Straight 5mm blade width.

ELX5ST

Periostomes

Designed to cut ligament fibers during the extraction procedure. Excellent access is achieved because the design incorporates the optimum blade angles. The DuraLite® ROUND handle has a large 3/8" (9.5mm) diameter and is the lightest handle available. This makes manipulation of the instrument unsurpassed for control and comfort. Features: Long, Slender, Ultra-sharp blades with two cutting edges.



Periostome #1 Anterior with a straight blade and socket expander.

REPTOME1



Periostome #2 Contra-Angle for posteriors.

REPTOME2