

Used to contour and remove alveolar bone during extraction procedures. Nordent Rongeurs are hardened and tempered to maintain their sharpness and have a satin finish to eliminate glare.



**Blumenthal #2**  
Standard profile blades cut around the entire periphery and are set at a 30° angle. (6-1/4" / 160 mm).

R2



**Mini-Friedman**

Narrow profile blades cut around the entire periphery and are set at a 30° angle. (5" / 125mm).

R1A



**Friedman #1**

Narrow profile blades cut around the entire periphery and are set at a 30° angle. (5-1/2" / 140 mm).


R1



**Niro Side Cutting #3**

Single straight cutting surface blade. This rongeur has a slight "S-Curve" angulation for easier access. (5-1/2" / 140 mm).

R3



**Cleveland #4**

Standard profile blades are set at a 30° angle. (5-1/2" / 140 mm).

R4