

# Root Tip Picks

Root Picks are used to retrieve loose root fragments after extraction. Nordent uses a special heat treating process that assures a sharp long lasting point and shank strength. The handle is all stainless steel and can be sterilized by any method.

Double ended Root Tip Picks and Root Teasers simplify tray set ups for more efficient procedures. Nordent offers these instruments in large diameter DuraLite® ColorRings™ handle or in the Standard handle.



**Heidbrink #1**  
Straight.

**Heidbrink #2**  
Left.

**Heidbrink #3**  
Right

**Heidbrink #2-3**  
Double ended root pick with pointed tips.

**Davis #11**  
Root tip teaser with a rounded tip.

**Davis #11P** Root tip teaser with a pointed tip and serrated blade.

ER1 shown  
ER1SR (serrated)

ER2 shown  
ER2SR (serrated)

ER3 shown  
ER3SR (serrated)

Handle Selection  
 CEER2-3 shown  
 ER2-3

Handle Selection  
 CEER11 shown  
 ER11

Handle Selection  
 CEER11P shown  
 ER11P