

Scissors

(shown 50% actual size)

Iris



Iris #301

Straight blades
(4 1/2" / 115 mm).

S301



Iris #302

Curved blades
(4 1/2" / 115 mm).

S302



Iris #303

Angular blades
(4 1/2" / 115 mm).

S303

Castroviejo



Castroviejo #311

Curved blades with spring action handle. Excellent for tissue and suture. (4" / 100 mm).

S311

Suture



Suture #304

Angular blades with suture hook.
(4 1/2" / 115 mm).

S304



Spencer Suture #306

Straight blades with suture hook.
(4 3/4" / 120 mm).

S306



Suture #322

Long handle, straight blades with suture hook.
(6 1/4" / 160 mm).

S322

LaGrange



LaGrange #314

A unique curvature allows easy access to the tightest areas. One blade is serrated to hold tissue when cutting. (4 1/2" / 115 mm).

S314

(shown 50% actual size)

Goldman-Fox

Sullivan



Goldman Fox #313
Curved blades with one blade serrated to hold tissue when cutting. (5" / 130 mm).

S313



Goldman Fox #321 A unique handle curvature allows easy access to the tightest areas. One blade is serrated to hold tissue when cutting. (5" / 130 mm).

S321




Sullivan #316 Extra long curved blades. One blade is serrated to hold tissue when cutting. (5 1/2" / 140 mm).

S316


Kelly

Dean



Kelly #317
Long handle with straight blades. (6 1/4" / 160 mm).

S317




Kelly #318
Long handle with curved blades. (6 1/4" / 160 mm).

S318



Kelly #319
Long handle with angular blades. (6 1/4" / 160 mm).

S319



Dean Gum #320 Long curved handle with straight blades that are angled upward. One blade is serrated to hold tissue when cutting. (6 3/4" / 170 mm)

S320